


Dear Customer,

From our experience, VS range ventilators can be used in ventilation in the following configurations provided that patients can use the ventilator correctly, that it is fixed properly to keep from falling and that the alarms have been set accurately.

However, it shall be used within the bounds of common sense including activities of daily living, and it shall be powered by the prescribed external battery or internal battery.

Allowed configurations:

1. during transport in hospital wards on a stretcher or in a wheelchair,
2. during transport with a helper, carrying VS in a bag,
3. during transport in a wheelchair,
4. during transport in a private car to go to hospital,
5. on flights with permission from airlines.



2008/08/29

Britta JAFFRE

*Saime – ResMed Group  
Quality & Regulatory Director*

Saime SAS, société du groupe ResMed

Z.I - 25 rue de l'Etain - 77176 Savigny-le-Temple - France - Tél.: +33 (0)1 64 19 11 11 - Fax: +33 (0)1 64 41 81 30

Société par actions simplifiée au capital de 7 890 000 € - Siret Melun 451 683 536 00028 - APE/NAF 741J - N° TVA CEE FR02451683536